



New Upgrade ! Venticap 502 Filtration Box

Bigger ! Stronger Fan ! More efficient ! with the same Price



Venticap H502



New Venticap, New Name :
Chemtrap H402 SE





BENEFIT

- **Safe!** – Filtration efficiency compliance to NFX 15-211, Class 2 standard
- **Ventilation** – provide more than 10 air change per hour, according to EN 14770-01
- **Easy, Universal** – Easy installation and fits to all type of flammable & acid storage cabinet
- **Protection of User' Health** – No inhalation of the toxic chemicals when opening the cabinet
- **24 hours Air Cleaning Function!** – 24/7 cleaning of air of the room at a rate of 75 m³/hr
- **Cost Saving!** – No connection to HVAC system, save energy bills and reducing greenhouse gases.

4 mains upgrade :

- Stronger Air Flow : 66 m³/hr
- Quieter with only 50 dbA
- Easy Filter Replacement without metal screw
- With option of MF2 filter and dual F2 filters

DIMENSIONS



Product Name	Chemtrap H402 Midcap	Chemtrap H402 SE (New Venticap)	Venticap 502
External Width(mm)	482	450	570
External Depth(mm)	451	450	390
External Height(mm)	315	335	270
Air Flow (m ³ /h)	60	66*	62*
Safety Standards	Filtration performance tested according to the AFNOR NFX 15-211:2009 STANDARD: FRANCE EN1822:1998-CE Marking	Filtration performance tested according to the AFNOR NFX 15-211:2009 STANDARD: FRANCE EN1822:1998-CE Marking	Filtration performance tested according to the AFNOR NFX 15-211:2009 STANDARD: FRANCE EN1822:1998-CE Marking
Voltage/Frequency (V/Hz)	100-240V/50-60Hz	100-240V/50-60Hz	230-240V/according to the country
Power Consumption	10W/0.4A	10W/0.4A	65W/0.8A
DC Output Voltage	24VDC	24VDC	N/A
Noise	50dbA	50dbA*	57dbA*
Filter Replacement	Lift the filter housing with the hydraulic rods	Remove the filter housing	Unscrew the screws needed