


ALARA STORAGE CABINET



PPE (PERSONAL PROTECTIVE EQUIPMENT) STORAGE CABINET

Under the Occupational Safety and Health Act (OSHA), employers are responsible for providing safe and healthful workplaces for their employees. Personal Protective Equipment (PPE) are required in every workplace to improve safety and health of their employees. Personal Protective Equipment (PPE) are to be stored properly and easily accessible by the employees. PPE storage cabinet is essential to provide large space to store medical, laboratory and industry used Personal Protective Equipment (PPE). ALARA's PPE storage cabinet made of high quality acid corrosion powder coated steel. Designed for durability and user-friendliness.

PPE Storage are used to store PPE for (1) Eye and Face Protection, (2) Foot Protection (3) Head Protection, (4) Protective Clothing and Gloves (5) Knee and Elbow Protection (6) Respiratory Protection, (7) Hearing Protection and (8) Fall Protection.




PPE STORAGE CABINET

Brand	ALARA, Malaysia
Model	T-PC1020-MW
Dimensions	920W x 500D x 1800H mm

- Left Compartments with 3 adjustable Shelves
- Right Compartment with 1 Hanging
- Comes with 2 Swing Doors with Wire Mesh
- Blue Colour
- 0.8mm thickness epoxy coated steel

MICROSCOPE STORAGE CABINET

ALARA Microscope Storage Cabinet is made with chemical resistant powder coated steel. Cabinet Model T-MC1020-MG stores up to 12-16 microscopes. All items are kept neat and organised on the 3 adjustable shelves, large space and protected by 2 Swing Doors with lock. Microscope storage cabinet is build in with temperature controller to achieve the required humidity in the storage chamber. This cabinet is an effective way to store microscopes and prevent fungus growth on the microscope lenses. Colour options are available and customisable.



MICROSCOPE STORAGE CABINET

Brand	ALARA, Malaysia
Model	T-MC1020-MG
Dimensions	900L x 450D x 1960H mm

- 2 Swing Doors with Insert glass & lock
- Comes with 3 adjustable shelves
- Build in temperature controller complete with sensor & interfacing circuit





LABORATORY SAMPLE STORAGE CABINET

The ALARA's laboratory sample storage cabinet is constructed with thick high quality chemical resistant powder coated steel. ALARA's laboratory sample storage cabinet has a large 900 mm width x 450 mm depth with height of 1750 mm, and comes with 3 adjustable shelves. Customer could choose various type of door design : from open type storage cabinet, 2 swing metal doors, 4 swing metal doors to full glass doors. This will allow customer to mix and match the sample cabinets in the laboratory to provide a modern, clean and aesthetic environment.

OPEN TYPE STORAGE CABINET

Brand	ALARA, Malaysia
Model	T-AS1020
Dimensions	900W x 450D x 1750H mm

- Comes with 3 adjustable shelf



STORAGE CABINET WITH 2 SWING METAL DOOR



Brand	ALARA, Malaysia
Model	T-AS1020-M
Dimensions	900W x 450D x 1750H mm

- Comes with 3 adjustable shelf

STORAGE CABINET WITH 2 SWING INSERT GLASS DOOR

Brand	ALARA, Malaysia
Model	T-AS1020-MG
Dimensions	900W x 450D x 1750H mm

- Comes with 3 adjustable shelf
- Doors comes with 4 insert glass for visible purpose



STORAGE CABINET WITH 4 SWING METAL DOOR

Brand	ALARA, Malaysia
Model	T-AS1220-M
Dimensions	900W x 450D x 1750H mm

- 2 Tiers with 4 Swing Metal Door
- Each tier comes with 2 adjustable shelf (Total 4 adjustable shelf)



ALARA STORAGE CABINET



STORAGE CABINET WTH 2 SWING INSERT GLASS DOOR

Brand	ALARA, Malaysia
Model	T-AS1220-MG
Dimensions	900W x 450D x 1750H mm

- 2 Tiers with 4 Swing insert Glass Door for visible purpose"
- Each tier comes with 2 adjustable shelf (Total 4 adjustable shelf)



STORAGE CABINET WTH 2 SLIDING GLASS DOOR



Brand	ALARA, Malaysia
Model	T-AS1220-LG
Dimensions	900W x 450D x 1750H mm

- 2 Tiers with 2 Sliding Glass Door for visible purpose
- Each tier comes with 2 adjustable shelf (Total 4 adjustable shelf)

STORAGE CABINET WTH 2 SWING INSERT GLASS DOOR AND 2 SWING SOLID DOOR

Brand	ALARA, Malaysia
Model	T-AS1220-MGM
Dimensions	900W x 450D x 1750H mm

- 1 Tier with 2 Swing Insert Glass Door for visible purpose
- 1 Tier with 2 Swing Solid Door
- Each tier comes with 2 adjustable shelf (Total 4 adjustable shelf)



STORAGE CABINET WITH 1 SLIDING GLASS DOOR AND 1 SLIDING SOLID DOOR


Brand	ALARA, Malaysia
Model	AS 1220 LGLM
Dimensions	900W x 450D x 1750H mm

- 1 Tier with 1 Sliding Glass Door for visible purpose
- 1 Tier with 1 Sliding Solid Door
- Each tier comes with 2 adjustable shelf (Total 4 adjustable shelf)

STORAGE CABINET WTH 2 SWING INSERT GLASS DOOR AND 2 SWING SOLID DOOR

Brand	ALARA, Malaysia
Model	AS 1220 LGMG
Dimensions	900W x 450D x 1750H mm

- 1 Tier with 1 Sliding Glass Door for visible purpose
- 1 Tier with 1 Sliding Metal Insert Glass Door
- Each tier comes with 2 adjustable shelf (Total 4 adjustable shelf)



TEAM MEDICAL & SCIENTIFIC SDN BHD

Unit D2-5-G (Ground Floor), Block D2, Pusat Perdagangan Dana 1, Jalan PJU 1A/46, 47301 Petaling Jaya, Selangor, Malaysia.

T. +603-7832 0736 | F. +603-7831 3210 | E. account@tms-lab.com | W. www.tms-lab.com